

Stamina Antibody Dilution Buffer, Cat.# R2004

- **2x more sensitive:**
 - comparing to using skim milk
- **Convenient:**
 - use at room temperature
- **Save antibody:**
 - re-use the antibody 3-6 times

Instruction

Storage: 4°C

1. Add your primary or secondary antibody directly in this buffer at a desired dilution. Submerge the membrane in the buffer. It is not necessary to block the membrane with milk or BSA prior to this step.
2. Incubate the blot in the buffer for 1-24 hours at **room temperature**. Longer incubation time will linearly increase the final signal strength.
3. Recover the primary antibody solution and store at 4°C
4. Proceed the downstream steps as described in product R1100 protocol.